



Ensuring Nucleic Acid Quality for Reliable Research Outcomes

Demonstrating Versatile Performance of the Implen NanoPhotometer® C40

The success of subsequent procedures like PCR, gene expression analysis, and metagenomic sequencing is directly impacted by accurate nucleic acid quantification, which is a crucial stage in molecular biology and environmental research. This paper reports two practical applications of the Implen NanoPhotometer® C40: siRNA-based gene silencing in breast cancer research and environmental microbiology, where DNA was tracked during petroleum biodegradation. The NanoPhotometer® C40 ensured data quality and reproducibility in both situations by enabling quick, accurate DNA and RNA quantification with small sample volumes. The instrument is a vital tool in a variety of scientific contexts due to its sturdy design, full-spectrum detection, and contamination alerts.

Keywords or phrases:

NanoPhotometer C40, nucleic acid quantification, UV-Vis Spectrophotometry, DNA purity, RNA quality control, environmental bioremediation, siRNA gene silencing, qRT-PCR validation, microbial community profiling, Implen spectrophotometer

INTRODUCTION

In workflows ranging from therapeutic gene expression studies to microbial community profiling, accurately and confidently quantifying nucleic acids is essential. Whole experiments may be jeopardised by inconsistent measurements or contaminants that are not yet detected. Even though they are accurate, traditional UV-Vis instruments frequently need larger sample volumes, cuvettes, and labour-intensive maintenance. In order to overcome these difficulties, the Implen NanoPhotometer® C40 provides full-spectrum, ultra-low volume quantification without the need for cuvettes or recalibration, as well as immediate contamination alerts.

This study demonstrates the use of the NanoPhotometer® C40 in two research settings:

(1) determining the integrity of DNA during hydrocarbon bioremediation, and (2) analysing the quality of RNA in a model of breast cancer cells for RNA interference investigations.

INSTRUMENT OVERVIEW

Protein and nucleic acid concentrations can be determined using the Implen NanoPhotometer® C40, a UV-Vis spectrophotometer, from volumes as small as 0.3 µL. It allows for high dynamic range analysis without sample dilution by supporting pathlengths of 0.67 mm and 0.07 mm. With integrated performance verification and no need for recalibration, the system offers automatic contamination detection and purity evaluation using absorbance ratios (A260/280 and A260/230).

Case Study 1:

DNA Quantification in Bioremediation

The ability of enriched microbial communities to break down crude oil in marine sediment was examined in the study "Bioremediation of Petroleum Hydrocarbons Using Microbial Consortia." To improve the microbial consortia's capacity to degrade hydrocarbons, they were grown in oil-supplemented conditions during several enrichment cycles.

The MO BIO PowerSoil® DNA Isolation Kit was used to extract genomic DNA from culture samples after each enrichment. Before moving on to microbial profiling procedures like PCR-DGGE (Denaturing Gradient Gel Electrophoresis) and Illumina sequencing, the NanoPhotometer® C40 was used to evaluate the quantity and quality of DNA.

Experimental

DNA concentration was measured to guarantee uniformity among enriched consortia.

To verify sample purity and rule out the possibility of protein or salt contamination, the A260/280 and A260/230 ratios were employed.

Researchers were able to identify any extraction carryover

(such as leftover phenol or humic acids), which is typical in soil-based samples, thanks to real-time spectrum visualisation.

The C40's precise quantification made sure that only high-quality DNA moved on to sequencing, boosting the accuracy of analyses of microbial diversity and cross-enrichment generation comparisons. As a result, strong consortia dominated by *Alcanivorax*, *Marinobacter*, and other hydrocarbon-degrading taxa were identified.

Case Study 2:

RNA Quantification in Breast Cancer Research

In order to silence the gene CRIF1 in MCF-7 human breast cancer cells, researchers created a siRNA-loaded lipid nanoparticle (LNP) system, which was published in the *International Journal of Molecular Sciences*. The objective was to assess how CRIF1 knockdown affected p53, p21, and TIGAR-mediated tumor-suppressive pathways.

Following LNP transfection, TRIzol reagent was used to extract total RNA from both treated and control cells. Prior to reverse transcription and qRT-PCR, the extracted RNA was examined with the NanoPhotometer® C40.

Experimental

RNA concentrations were measured in order to standardise the inputs for cDNA synthesis across conditions. RNA purity appropriate for enzymatic reactions was confirmed by A260/280 ratios greater than 2.0.

Degradation or leftover solvents could be found using visual absorbance profiling, which is essential for maintaining the accuracy of subsequent gene expression analysis.

The researchers confidently carried out the qRT-PCR analysis, which demonstrated successful gene silencing and related downstream effects, because of the NanoPhotometer®'s accuracy.

RESULTS AND DISCUSSION

Because it provided quick, accurate nucleic acid quantification with small sample volumes and no consumables, the NanoPhotometer® C40 was essential to both investigations. The highest analytical standards were maintained in clinical and environmental research workflows thanks to its capacity to compute purity ratios and provide instant visual feedback.

It enabled the following in the bioremediation study:

- ▶ DNA input that is constant throughout DGGE and sequencing runs
- ▶ Dependable observation of the dynamics of microbial communities

It made it possible for the oncology study to:

- ▶ Normalisation of RNA input with confidence for qRT-PCR
- ▶ Reliable confirmation of the therapeutic effects of siRNA

CONCLUSION

In two different fields—microbial environmental monitoring and molecular oncology—the Implen NanoPhotometer® C40 showed great utility. It is a vital tool in contemporary labs because of its capacity to provide accurate nucleic acid quantification results with incredibly small sample volumes. The C40 guarantees that sample quality is never a bottleneck, whether it is for analysing RNA from cancer cells or sequencing bacterial DNA from sediments contaminated by oil.

REFERENCES

1. Sandia Primeia; Inoue, C.; Chien, M.-F. Potential of Biosurfactants' Production on Degrading Heavy Oil by Bacterial Consortia Obtained from Tsunami-Induced Oil-Spilled Beach Areas in Miyagi, Japan. *Journal of Marine Science and Engineering* 2020, 8 (8), 577–577. <https://doi.org/10.3390/jmse8080577>.
2. Piao, S.; Lee, I.; Kim, S.; Park, H.; Nagar, H.; Choi, S.-J.; Vu, G.-H.; Kim, M.; Lee, E.-O.; Jeon, B.-H.; Kim, D. W.; Seo, Y.; Kim, C.-S. CRIF1 siRNA-Encapsulated PLGA Nanoparticles Suppress Tumor Growth in MCF-7 Human Breast Cancer Cells. *International journal of molecular sciences* 2023, 24 (8), 7453. <https://doi.org/10.3390/ijms24087453>.